

E.C.G Electrodes

Liquid Gel Electrodes

Fast ECG signal transmission even if inadequate clean skin, give dependable quick monitoring with easy application and removal.

S-55 (Foam backing , Ø 55mm)

S-50 (Foam Backing , Ø 50mm)

M-55 (Micro porous Tape , Ø 55mm)

M-50 (Micro porous Tape , Ø 50mm)



E.C.G Electrodes

Solid Gel Electrodes

In this type Gel dry out is not an option, it also can be pressed in the centre without the risk of gel spreading onto the adhesive materials. Due to gel character, it is very good choice for long term monitoring and reduction of motion artifacts.

MD-50 (Micro porous Tape , Ø 50mm)

CD-50 (Soft Cloth , Ø 50mm)

SD-50 (Foam Backing , Ø 50mm)

SD-45 (Foam Backing , Dimension 45 x 50mm)

SD-40 (Foam Backing , Dimension 40 x 45mm)

SD-36 (Foam Backing , Dimension 36 x 50mm)

SD-30 (Foam Backing , Dimension 30 x 35mm)

SD-26 (Foam Backing , Dimension 26 x 49mm)

T-tracet (Foam Backing , Dimension 36 x 55mm)

FSD-36 (Enhanced Foam Backing for Stress Test , Dimension 36 x 50mm)



Special Electrodes

Electrodes with banana plug, Neonatal Electrodes with pre-attached cables and Radio-translucent electrodes

