

# Spinal Needle

Key/slot arrangement of stylet and cannula hubs facilitates proper needle bevel orientation Unique translucent window hub featuring contact clarity that helps allow visualization of CSF Fitted stylet reduces tissue coring.

## Standard Length

18G x 3 1/2" 23G x 3 1/2"  
19G x 3 1/2" 25G x 3 1/2"  
20G x 3 1/2" 26G x 3 1/2"  
21G x 3 1/2" 27G x 3 1/2"  
22G x 3 1/2"

